

**Works Better.  
Lasts Longer.  
Costs Less.**

**DEFEND<sup>®</sup>**  
by Mydent International



[www.defend.com](http://www.defend.com)





## The Next LEVEL In Face Mask Protection...

See Page 3 & 4



**PLACE YOUR ORDER THROUGH YOUR DEALER!**  
**[www.Defend.com/marketplace](http://www.Defend.com/marketplace)**

Mydent International is excited to launch its new DEFEND Marketplace! For a limited time, healthcare professionals who order through this website will receive a \$50 Restaurant.com eGift Certificate with orders over \$150.

- Order through your dealer of choice
- View Defend online promotions
- Quick order processing (usually ships next business day)
- View prior orders made on website
- Save your "favorite" products
- User friendly

**Use Coupon Code FREE50**

For your Restaurant.com eGift card.





**Works Better.  
Lasts Longer.  
Costs Less.**



## IT'S ALL HERE.

DEFEND is a complete line of infection control products, disposables, preventatives, impression material systems, rotary instruments, and restoratives. Backed by superior service, outstanding customer support and a global network of distribution partners, DEFEND delivers convenience, satisfaction and the highest quality products available today.

## SUPERIOR VALUE.

Unbeatable prices meet unmatched quality. By choosing the DEFEND brand, customers enjoy the exceptional value of products that **"Work Better. Last Longer. Cost Less."**

## EXPERIENCE MEETS INNOVATION.

With over thirty years of expert production and support, Mydent boasts unmatched knowledge of the healthcare supply industry. We recognize the continual need and unrelenting demand for aggressive innovation. From product breakthroughs to customer rewards programs, Mydent continuously strives to deliver the best products and services available.

## HERE. THERE. EVERYWHERE.

From the Americas to the Far East, Mydent has built a global presence through diligence, superior products and unsurpassed support. With key production and shipping facilities around the world and a growing network of importers and distributors, DEFEND products are moving faster than ever.

## ASK FOR FREE SAMPLES.

[www.defend.com/freesamples](http://www.defend.com/freesamples)

Mydent on Facebook:  
<http://www.facebook.com/MydentDefend>

Mydent on Twitter:  
[https://twitter.com/mydent\\_defend](https://twitter.com/mydent_defend)

Mydent on LinkedIn:  
<https://www.linkedin.com/company/mydent-international>

Mydent on Google Plus:  
<https://plus.google.com/+DefendMydent>

Mydent on Youtube:  
<https://goo.gl/h7Q29o>

NEW! Mydent on Instagram:  
<https://www.instagram.com/mydentinternational>

The information contained herein is provided for informational purposes only, without any representation or warranty of accuracy or completeness of information or other warranty of any kind. In particular, we can make no warranties, either expressed or implied, that the products in this catalog will fully protect the user from exposure to blood or body fluids or risk of contracting infectious diseases. Additionally, we may make improvements or changes in our products at any time, without notice. This catalog may not be copied in whole or in part without the express written permission of Mydent International, Inc. ©2017



## TABLE OF CONTENTS

FACE MASKS & SHIELDS	<b>3-4</b>
GLOVES	<b>5</b>
DISINFECTANTS & CLEANERS	<b>6</b>
EVACUATION PRODUCTS	<b>7</b>
STERILIZATION PRODUCTS	<b>8</b>
IMPRESSION MATERIALS	<b>9-10</b>
MIXING TIPS	<b>11</b>
PREVENTATIVES	<b>12</b>
DISPOSABLES	<b>13-14</b>
BARRIER PRODUCTS	<b>15</b>
INSTRUMENTS	<b>16</b>
RESTORATIVE	<b>16</b>
APPAREL & EYEWEAR	<b>17-18</b>



## Reward!

**Earn Reward Points for all DEFEND Products Purchased!**  
**Redeeming Gifts is Easy and Our Selection is Extraordinary.**

**For all NEW members sign up and receive:**

☒ **DOUBLE POINTS**  
and

☒ **1,000 FREE POINTS**

**Go to [defend.com/rewards-program](http://defend.com/rewards-program) to cash in on your earned DEFEND product points!**


All DEFEND masks are FREE from:



## Dual Fit Ear-Loop Face Masks (Pleated)

*Dual Fit chin contour to reduce gaps for maximum protection!*

**NEW!**



**Features**  
Pliable  
Chin Band!

### Level 1

#### Features & Benefits:

- 4 fold pleat
- Smooth cellulose inner layer for sensitive skin
- Soft, rounded ear loops
- Form-fitting pliable nose & chin band

50 per box  
MK-7100 Blue  
MK-7110 Pink  
MK-7120 White



**Replaces:**

#MK-1410, #MK-2000, #MK-2001, #MK-2007, #MK-2008 & #MK-2010



### Level 1

#### "Diffuser" Anti-Fog

#### Features & Benefits:

- Anti-fog, moisture absorbent nosepiece
- Soft, rounded ear loops
- Soft, 3-ply pleated material for ultimate comfort
- Form-fitting pliable nose & chin band

50 per box  
#MK-1046 Blue  
#MK-1056 Pink  
#MK-1076 Blue (Tie-On)




Close-up of Anti-Fog, Moisture Absorbent Nosepiece Strip

**1 LEVEL**

ASTM 1

Designed for procedures with low amounts of blood, fluid, spray and/or aerosol exposure.

**NEW!**



**Features**  
Pliable  
Chin Band!

### Level 2

#### Features & Benefits:

- 4 fold pleat
- Smooth cellulose inner layer for sensitive skin
- Soft, rounded ear loops
- Form-fitting pliable nose & chin band

50 per box  
MK-7200 Blue  
MK-7210 Pink  
MK-7220 White



**Replaces:**


#MK-4046

**2 LEVEL**

ASTM 2

Ideal for procedures where moderate to light amounts of fluid, spray and/or aerosols are produced.

**NEW!**



**Features**  
Pliable  
Chin Band!

### Level 3

#### Features & Benefits:

- 4 fold pleat • 4 Ply
- Smooth cellulose inner layer for sensitive skin
- Fluid resistant Outer Layer
- Soft, rounded ear loops
- Form-fitting pliable nose & chin band

50 per box  
MK-7300 Blue  
MK-7310 Pink  
MK-7320 White



**3 LEVEL**

ASTM 3

Designed for procedures with moderate to heavy amounts of blood, fluid, spray and/or aerosol exposure.



### "Breathe E-Z"

#### Dual Fit Ear-Loop Face Mask

- Economical
- Easy breathing
- Soft, rounded ear loops
- Soft inner lining for comfort
- Form-fitting pliable nose & chin band

50 per box  
MK-1246 Blue  
MK-1256 Pink  
MK-1266 Lavender

#### Face Mask Material Requirements by Performance Level

ASTM F2100-11 STANDARDS	ASTM Level 1	ASTM Level 2	ASTM Level 3
FLUID RESISTANCE, mmHG	80	120	160
BFE	≥95%	≥98%	≥98%
PFE, @0.1 micron	≥95%	≥98%	≥98%
DELTA P, mm H <sub>2</sub> O/cm <sup>2</sup>	<4.0	<5.0	<5.0
FLAME SPREAD	Class 1	Class 1	Class 1



Pediatric patients *and* hygienists will love our New Krazy Cone Face Masks featuring four delightful character designs.

### Features:

- Non-Glare Pliable Nosepiece
- Cool, Comfortable Breathing
- Fluid Resistant
- Soft Inner Layer for Facial Comfort
- Fiberglass & Latex Free
- 99%+ BFE & 99%+ PFE



50 per box  
 #MK-3100 Dog  
 #MK-3200 Cat  
 #MK-3300 Rabbit  
 #MK-3400 Pig  
 40 per box  
 #MK-3500 Assorted



### Premium Molded Mask

These fluid-resistant masks offer a superior BFE and PFE of greater than 99%. Available in both regular size and a smaller size in pink, each offering a non-glare, adjustable nosepiece for a comfortable fit.

50 per box  
 #MK-1006 Blue  
 #MK-1016 Pink (10% smaller)



### Premium Ear-Loop Molded Mask

Innovative ear-loop design cone mask loops around the ears instead of the back of the head, eliminating hair tangles. Designed for easy breathing, are highly pliable, fluid resistant, and features a non-glare nosepiece. Offers superior protection and a BFE/PFE of greater than 99%.

50 per box  
 #MK-1036 Blue  
 #MK-1066 Pink (10% smaller)



### "Breathe-EZ" Valueline Molded Mask

These comfortable molded masks are the easiest to breathe through. Generally used for examinations and any non-splatter, dustless, or fluidless procedures. Offers superior protection and a BFE/PFE of greater than 99%.

50 per box  
 #MK-1106 Blue



### Economy Molded Mask

This lightweight mask offers a pliable nosepiece with a single band for a cool and comfortable fit. Offers good protection at a BFE/PFE of greater than 99%.

50/box 100/box  
 Blue #MK-1007 #MK-1000  
 Pink #MK-1008 #MK-1001



### Form-Fit Mask

The unique contour design without a metal nosepiece prevents fogging of eyewear. These masks are ergonomically designed, fluid resistant and offer superior protection and a BFE/PFE of greater than 99%.

50 per box  
 #MK-1206 Blue  
 #MK-1216 Pink



### N95 Particulate Respirator

A NIOSH approved N95 Particulate Respirator that meets guidelines for TB exposure control. Indicated for use when treating patients with airborne diseases such as TB or Influenza.

20 per box  
 #MK-N95



### Level 3 Pleated Ear-Loop Mask w/Shield

- 4 Fold Pleat
- Smooth cellulose inner layer for sensitive skin
- Dual Fit chin contour to reduce gaps for maximum protection
- Soft, rounded ear loops
- Form-fitting pliable nose & chin band
- Fiberglass/latex free
- Optically Clear Shield

25 per box  
 #MK-7400 Blue



### Face Shields

This face shield is made to fit over glasses and face masks. The adjustable velcro headband provides a comfortable fit. The shield protects against splashing. It is anti-fog treated and latex and fiberglass free. Available in full and 1/2 shield sizes.

25 per box  
 #MK-3000 Full Shield  
 #MK-3001 Half Shield





## Powder-Free Micro-Textured Latex Exam Gloves

Our Powder-Free Latex Exam Gloves are thin for extra tactile sensitivity and have a micro-textured finish. The gloves are polymer coated for easy on/off comfort.

100 per box

#LG-2102 X-Small (6)

#LG-2103 Small (6 ½)

#LG-2104 Medium (7)

#LG-2105 Large (8 ½)

#LG-2106 X-Large (9)



## Blackjack Powder-Free Textured Latex Gloves

Blackjack Powder-Free Textured Latex Gloves are polymer coated for easy on/off. The textured finish creates an enhanced grip.

100 per box

#LG-8002 X-Small (6)

#LG-8003 Small (6 ½)

#LG-8004 Medium (7 ½)

#LG-8005 Large (8 ½)

90 per box

#LG-8006 X-Large (9)



## Powder-Free Vinyl Exam Gloves

Our Powder-Free Vinyl Exam Gloves are made of high-quality stretch vinyl, are tear resistant and have an extra soft feel. Vinyl gloves are a great alternative for individuals who may be sensitive to latex.

100 per box

#VG-4003 Small (6 ½)

#VG-4004 Medium (7 ½)

#VG-4005 Large (8 ½)

#VG-4006 X-Large (9)



## Powder-Free Nitrile Textured Exam Gloves

Powder-Free Nitrile Textured Exam Gloves are highly puncture resistant, have great tactile strength and dexterity and are a great alternative for individuals sensitive to latex materials.

100 per box

#NG-2002 X-Small (6)

#NG-2003 Small (6 ½)

#NG-2004 Medium (7 ½)

#NG-2005 Large (8 ½)

#NG-2006 X-Large (9)





## DEFEND + PLUS Disinfecting/Deodorizing/Cleaning Wipes

These disinfecting wipes contain a stable, low pH formulated disinfectant and deodorant for use on hard, non-porous surfaces in hospitals, emergency medical settings, laboratories, medical offices, dental offices – anywhere control of cross contamination is required. Containing a dual chain quaternary ammonium/ alcohol formula, the product is proven effective against the following pathogens in **2 minutes\*** (see below) at 20°C/60°F, on hard, non-porous surfaces when used as directed.

#SO-9000 - 160 per canister - 6"x 6" sheets  
#SO-9001 - 60 per canister - 10"x 10" sheets

- **New expanded 2 minute kill time list\***
- **Disinfects, cleans & deodorizes**
- **Kills 99% of germs on contact**
- **EPA registered hospital grade disinfectant**
- **Helps prevent cross contamination**
- **Convenient and ready to use with no mess**

- ★ • *Escherichia coli* 0157:H7
  - Human Hepatitis B virus
  - Human Hepatitis C virus
  - Herpes Simplex virus 2 (genital herpes virus)
  - HIV-1 (AIDS virus) - 1 minute kill [See label for complete information]
  - Influenza A (A2 Japan - H1N1)
  - Methicillin Resistant *Staphylococcus aureus* (MRSA)
  - *Mycobacterium tuberculosis*
  - *Pseudomonas aeruginosa*
  - *Salmonella enterica*
  - *Staphylococcus aureus*
  - Vancomycin Resistant *Enterococcus faecalis* (VRE)
- 5 minute kill time:
- *Trichophyton mentagrophytes* (Athlete's foot fungus)



## Naturally Hydrating Lotion Soap

This gentle and non-drying soap is specially formulated to care for and protect frequently washed hands. It is specifically blended with Aloe Vera Gel and "BioStructure™" - a blend of 7 natural ingredients devised to moisturize and preserve. Has a pleasant baby powder scent.

1 Gallon Bottle #SO-9300      16 Oz. Bottle #SO-9310



## Autoclave Cleaner

Concentrate: 4 oz makes 68 ounces solution for steam autoclave cleaner, enhanced cleaning power, contains no phosphates. Cleans grease, grime, and particles from all sterilizers.

32 Ounce Bottle  
#AC-8000

## Ultrasonic Solutions



### General Purpose Cleaner (#1)

Non-ammoniated, 10:1 concentrated formula effectively cleans dried blood, tissue, pumice, food, rouge, Tripoli & buffing compounds from the surface of any appliance while providing a rust inhibitor.

1 Gallon Bottle - #SO-9400



### Tartar & Stain Remover (#4)

Use this potent stain remover on dentures, bridges, crowns and orthodontic appliances. Ready-to-use formula quickly breaks up and removes calculus and tartar while eliminating stains caused by tobacco and food.

1 Gallon Bottle - #SO-9600



### Temporary Cement Remover (#6)

This ready to use solution requires no mixing. The strong, fast-acting formula removes cement and stone quickly and completely from instruments, trays and other appliances.

1 Gallon Bottle - #SO-9800



**NOW WITH  
IMPROVED  
MULTI-ENZYME  
BLEND!**

## "SURG" Evacuation System Cleaner

Used daily, this superior evacuation system cleaner will keep lines running free and clear at maximum suction. Clinically proven to clean lines and discourage build-up. The non-foaming formula has a neutral pH for safe use with amalgam separators and maintains the ability to dissolve organic fluids and debris. The pleasant peppermint scent help to eliminate odors. 128 treatments per gallon.

1 gallon bottle  
(pump not included)  
#SO-9100



## "SURG" Dispenser Pump

In order to limit waste that goes into our landfills! Please request the SRG dispenser pump if you need one! Contact us today to receive the pump at no charge.

#SO-PUMP

## Trap-EZ Disposable Traps

DEFEND Trap-EZ disposable filters are designed to replace the trap presently in your evacuation unit or cuspidor. They reduce the possibility of infection and contamination, while improving the performance of your evacuation system.

144 per box



### #DT-5501

Fits A-dec models with black plastic canisters, plus newer models by Ampco, Belmont, DCI, Forest, Marus, MDT, & Schein - diameter 2 1/8".


### #DT-5502

Fits A-dec models with metal canisters - diameter 2 15/16".

### #DT-5503

Fits Ampco, Belmont, Biotec, Chayes, Dabi-Atlante, Dansereau, DCI, Den-Tal-EZ, Dentech, Forest, Knight, Pelton & Crane, Proma and others - diameter 1 7/8".

**\*Universal Diameter**



### #DT-6200C\*

Fits Siemens, Planmeca and many cuspidors with small diameter openings.

### #DT-HVE\*

Fits as a "cap" on all standard high volume evacuator tips.

### #DT-5508S\*

Fits all Saliva Ejectors.

### #DT-5509

Fits DCI/Belmont Reality models and System 6000 models - diameter 1 1/2"

### #DT-6100C\*

Fits most cuspidor bowls including: Adec, Ampco, Belmont, Dabi-Atlante, Dansereau, DCI, Den-Tal-EZ, Forest, Knight, Marus, MDT, Proma, Ritter & Schein.



**Clinically Proven!**

## Ultrasonic Enzymatic Tablets

Easy-to-use tablets dissolve in a gallon of water in under a minute to create a powerful, effective cleaning solution. Greatly assists in the removal of blood, bone, tissue, food particles and other soils that may be resistant to chemical detergents. Contains fast-acting protease enzymes and has a pleasant scent. Also use as an instrument presoak or evacuation system cleaner.

- Safe for amalgam separators
- Biodegradable, neutral pH • Makes 32 gallons
- Non-chlorine • Non-corrosive

64 tablets per box  
#UT-1000



## Standard Saliva Ejectors

These pliable, shape-maintaining saliva ejectors offer a comfortable tip that minimizes tissue impingement. Fits all standard saliva ejector hose ends.

100 per bag

#SE-7000 - White w/white tip

#SE-7001 - Clear w/blue tip

#SE-7002 - Assorted Colors



## Vented Oral Evacuation Tips

Convenient, single-use disposable tips contain smooth edges for maximum patient comfort. One end is vented and the other is non-vented. Fits all standard evacuator valves. Non sterile.

100 per bag


#VE-1121- 6" length



## Disposable Soft Touch Oral Evacuator Tips

- Cushioned, flexible and patient friendly
- Gentle on soft tissue and bone
- Reduces soft tissue aspiration
- Reduces noise

50 per bag  
#VE-2000



## Surgical Aspirator Tips

Made with smooth edges for maximum patient comfort and molded at a 30° angle - fits 11mm suction holders. Autoclavable up to 135°C or may be disposed after single use.

25 per bag

#ST-1020 - Blue 1/16" dia.

#ST-1021 - White 1/8" dia.


#ST-1023 - Green 1/4" dia.







- Thumb notch for easy opening of pouch
- Heavily sealed and reinforced bottom for added protection against instrument breaching
- Triple heat-sealed pouch to prevent possible penetration and tearing from sharp instruments
- Medical grade paper that meets or exceeds medical pouch standards
- Built in internal and external indicators
- Transparent film pouch for easy viewing of contents



## DEFENDLOC Sterilization Pouches with Built-in Internal Indicator



DEFEND+PLUS sterilization pouches are manufactured with lead-free, built-in dual internal and external indicators to ensure the correct sterilization temperature is reached in the autoclave chamber as well as inside the instrument compartment. Eliminates the need for internal indicator strips.

200 per box

#SP-0150 2.25" x 4"  
#SP-0250 2.25" x 2.75"  
#SP-0350 2" x 7.75"  
#SP-0550 2.75" x 10"  
#SP-1500 3.5" x 5.25"

#SP-2500 3.5" x 10"  
#SP-3500 5.25" x 10"  
#SP-4500 7.5" x 13"  
#SP-5500 9" x 15"  
#SP-8000 13" x 20"

1000 per box

#SP-2150 3.5" x 10"



## Sterilization Wrap (CSR)



Great for general purpose sterilization, provides a barrier against air and water-borne bacteria. Conforms easily around various shapes and sizes of cassettes, allowing excellent penetration of various methods of sterilization. Compatible with steam, ETO and gamma. Petroleum and memory-free material with a cloth-like feel.



1000/box #SW-1200 12" x 12"  
500/box #SW-2000 20" x 20"  
250/box #SW-3000 30" x 30"

1000/box #SW-1500 15" x 15"  
500/box #SW-2400 24" x 24"

## Sterilizer Indicator Tape

DEFEND Autoclave Tape is designed for use with autoclave or chemical sterilizers and will seal autoclave bags, tubing, pouches and CSR wraps. The tape will indicate when sterility temperature has been reached (275°F / 135°C). Made of extra tacky material, easily affixes to all types of paper and plastic.



60 yards per roll  
#AT-2001 1/2" inch  
#AT-2002 3/4" inch  
#AT-2003 1" inch



## Heat Sealer & Accessories

This ETL-certified and ultra durable heat sealer will seal nylon sterilization tubing within 1-2 seconds and then turns off automatically. Strong stainless steel, 8" heating elements resist breaking.



#HS-3000 110 Volt Unit  
#HS-3001 Heating Elements (3/pk)  
#HS-3002PK Upper Insulating Cloth (3/pk)  
#HS-3003PK Lower Insulating Cloth (3/pk)  
#HS-3004 Silicone Strip (opaque plastic)  
#HS-3005 Heat Sealer Repair Kit - Contains: 1 Heating Element, 1 Upper Cloth and 1 Silicone Strip

## Nylon Sterilization Tubing

Manufactured from transparent nylon for easy viewing of contents, the sterilization tubing will remain pliable after heat sterilization (up to 425°F/218°C). The tubing is for use in autoclaves, dry heat or chemical sterilizers and may be sealed with a heat sealer. The tubing comes in a 100 foot roll and has an optional easy-to-read sterilization indicator.



100 foot roll  
#AT-1018 2"  
#AT-1019 3"  
#AT-1020 4"  
#AT-1021 6"



## Super Hydrophilic Vinyl Polysiloxane Impression Material (VPS)



High performance (HP) cartridge system available in regular and fast set in all viscosities. Excellent tear strength to maintain impression integrity upon removal and precise impressions no matter what technique you use. Exceeds ISO 4380 standard which will help eliminate remakes.

4 cartridges per box

Viscosity	Regular Set	Fast Set
Light	#VP-9003	#VP-9008
Medium	#VP-9004	#VP-9009
Monophase	#VP-9005	#VP-9010
Heavy	#VP-9007	#VP-9011

## Super Hydrophilic Vinyl Polysiloxane Putty (VPS)

2 x 300ml jars

#VP-8006 - Regular Set

#VP-8008 - Fast Set



### VPS Impression Materials Setting Times

Item#		Viscosity	Color	Work Time	Set Time	Total Set Time
VP-9003 VP-9008	Light Body	Regular Set Fast Set	Yellow	2:30 1:00	1:00 1:00	3:30 2:00
VP-9004 VP-9009	Medium Body	Regular Set Fast Set	Buff	2:30 1:00	1:00 1:00	3:30 2:00
VP-9005 VP-9010	Monophase Body	Regular Set Fast Set	Purple	2:30 1:00	1:00 1:00	3:30 2:00
VP-9007 VP-9011	Heavy Body	Regular Set Fast Set	Rose	2:30 1:00	1:00 1:00	3:30 2:00
VP-8006 VP-8008	Putty	Regular Set Fast Set	Dk Blue	2:00 1:00	2:00 2:00	4:00 3:00



## Bite Registration Material



Mousse-like, fluffy vinyl polysiloxane material used for bite registration, copying transfers, and primary impression with bite trays. Thixotropic material has an excellent durometer of over 80 and flows easily, yet will not run or sag.

Contains: 2 cartridges (50 ml. each) plus 6 teal mixing tips

Unflavored

#BR-9001 Regular Set

#BR-9002 Fast Set

#BR-9003 Super Fast Set

Mint Chocolate Chip

#BR-9004 Fast Set

#BR-9005 Super Fast Set

"HARD" Bite – Unflavored

#BR-9010 Fast Set

### Bite Registration Materials Setting Times

Item#	Viscosity	Color	Work Time	Set Time	Total Set Time
BR-9001	Bite	Regular Set	2:00	1:00	3:00
BR-9002	Registration	Fast Set	:30	1:00	1:30
BR-9003	(Unflavored)	Super Fast Set	:30	:30	1:00
BR-9004	Bite	Fast Set	:30	1:00	1:30
BR-9005	Registration (Mint Choc Chip)	Super Fast Set	:30	:30	1:00
BR-9010	Hard Bite (unflavored)	Fast Set	:30	1:00	1:30
BR-9500	CADbite	Fast Set	:30	1:00	1:30
BR-9200	ClearBite	Fast Set	1:00	1:00	2:00



## Alginate Substitute

Polyvinylsiloxane "alginate alternative" impression material. Delivers exceptional accuracy and performance. Indicated for preliminary impressions, provisional crown and bridge impressions, study models, registration/opposing dentition impressions, impressions for orthodontic models, sports guards, bleaching trays and more.

- Work time 0:30 minutes
- Intra oral set 1:00 minute
- Total work/set time 1:30 minutes
- Durometer 80 (Shore A)
- Color: White



Contains: 2 cartridges (50 ml. ea.) plus 6 pink mixing tips  
#VP-9200



## Bite Trays

Single-use trays designed especially for the double-bite impression technique. Tray undercuts lock in impression material, facilitating accurate final impressions. Flexible, non-tear webbing and rigid frame design provide accurate and detailed impressions of both arches in a single procedure. Disposable dual-arch impression tray simultaneously captures master impression, counter impression and bite registration.

#BT-1001 Anterior 35/box #BT-1003 Sideless Post 50/box  
#BT-1002 Posterior 50/box #BT-1004 Quadrant 35/box



## Vinyl Tray Adhesive

Firm adhesion between impression trays and VPS impression materials. Available in 10ml w/ applicator brush.

#VP-8300 10ml




## ClearBite

Great visibility for voids or bubbles makes ClearBite an excellent choice as a matrix for quick and easy light cure temporary materials. Light cure through the clear material for dual cure or light cured temporary materials.

- Work time 1:00 minutes
- Intra oral set 1:00 minute
- Total work/set time 2:00 minutes
- Durometer 60-65 (Shore A)
- Clear & unflavored



Contains: 2 cartridges (50 ml. ea.) plus 6 pink mixing tips  
#BR-9200



## Impression Trays

Made of a strong, rigid material to assist in generating consistent impressions. Distortion free, preshaped, readily trimmed or flame softened for adjustments. Non-slip grips.

12 trays per bag

Perforated

#IT-7001 #1 large upper  
#IT-7002 #2 large lower  
#IT-7003 #3 medium upper  
#IT-7004 #4 medium lower  
#IT-7005 #5 small upper  
#IT-7006 #6 small lower  
#IT-7007 #7 quadrant ul/lr

#IT-7008 #8 quadrant ur/lr  
#IT-7009 #9 anterior upper  
#IT-7010 #10 anterior lower

Solid

#IT-7017 #7 quadrant ul/lr  
#IT-7018 #8 quadrant ur/lr  
#IT-7019 #9 anterior upper  
#IT-7020 #10 anterior lower





# Mixing Tips & Accessories



#VP-8204  
10:1 Gun

## Dispensing Gun

The dispensing gun is used to extrude impression materials and bite registration materials from cartridges. It consists of a release lever and a plunger, which are coordinated simultaneously to load and remove the cartridge. Works with 1:1 and 2:1 cartridges.

#VP-8202 For 1:1/2:1 HP Cartridge – 50ml

#VP-8204 For 4:1/10:1 Cartridge - 50ml  
(for use with Defend Crown & Bridge material)



## HP Style Mixing Tips



### HP Mixing Tips

High Performance Mixing Tips for all HP style 50ml impression material cartridges. Tips are color coded according to size.

48 per bag

#VP-8104 Teal 6.5mm  
#VP-8105 Yellow 4.2mm  
#VP-8106 Pink 5.4mm  
#VP-8107 Purple 7.5mm



## Crown & Bridge Mixing Tips

Tips are high performance, for temporary crowns, bridges and acrylics. Available in 10:1, 4:1 and 1:1 ratios.

25 per bag

#VP-8108 - 10:1 & 4:1 ratio

25 per bag

#VP-8111 - 1:1 ratio



## Temporary Cement Mixing Tips

Tips are high performance for temporary cement. Available in 1:1 ratio.

25 per bag

#VP-8140



## T-Mix Mixing Tips

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the New T-Mix tips is seamless: 100% compatible with the existing MIXPAC™ systems is guaranteed - as before, all original accessories fit.

48 per bag

#VP-8105T - Yellow  
#VP-8106T - Pink  
#VP-8104T - Teal



## Temporary Cement T-Mix Mixing Tips

Tips are high performance for temporary cement. Available in 1:1 ratio

25 per bag

#VP-8140T



## Ribbon Tips

For dispensing a ribbon-like layer of bite registration material on the occlusal surface. For use with mixing tip #VP-8104T (Teal).

24 per bag

#VP-8102



## Crown & Bridge T-Mix Tips


Tips are high performance, for temporary crowns, bridges and acrylics. Available in 10:1, 4:1 and 1:1 ratios

25 per bag

#VP-8108T - 10:1 & 4:1 ratio

25 per bag

#VP-8111T - 1:1 ratio



Recommended  
for Our Award  
Winning  
VPS & Bite!



## Core Material T-Mix Tips

Tips are high performance for core material. Also includes 25 clear intra-oral tips. Available in 1:1 ratio.

25 per bag

#VP-8150T



25 per bag

#VP-8165T



## Clear Intra-Oral Tips

Provide precise placement interproximally. Fits with universal mixing tips, or with T-Mix Mixing Tip #VP-8104T.

100 per bag

#VP-8002 Large

#VP-8101SM Small



## Yellow Intra-Oral Tips

Provide precise placement interproximally. Fits 4.2mm mixing tips and to be used with T-Mix Mixing Tip #VP-8105T.

100 per bag

#VP-8103



## DEFEND Desensitizer

### Universal Desensitizing Agent

- Quick & Easy Application
- Universal Compatibility
- Antimicrobial Agent

DEFEND Desensitizer is a superior desensitizing agent which can be placed under dental cements or temporary, provisional or final restorative materials to reduce post-operative sensitivity. It can also be used for desensitization of amalgam restorations, placed either conventionally or bonded.

DEFEND Desensitizer contains HEMA to aid bonding primers to penetrate etched dentin. The HEMA base in DEFEND Desensitizer helps seal tubules, giving immediate relief from pain due to hypersensitivity.



### APPLICATION IS QUICK AND EASY. DEFEND DESENSITIZER IS IDEAL FOR:

- Abrasion, cervical erosions and preps
- Applications to vital crown preparations prior to luting with cements
- At the "prep" appointment to desensitize during temporization

Most dentin resin bonding materials will benefit from DEFEND Desensitizer application. Its application reliably reduces post-op sensitivity by supporting the collagen framework for easier penetration of the adhesive, thus enhancing the dentin bond. DEFEND Desensitizer is an equivalent formula to market leader (for about 1/4 the cost per ml).

# DD-1000

10ml Bottle

HEMA, Glutaraldehyde

## Prophy Paste

Smooth, pliable and reduced-splatter formula contains fluoride and a unique blend of cleaning and polishing agents. Safely and effectively removes stains from tooth surfaces.



200 per box

	Coarse	Medium	Fine
Assorted*	#PP-1000	#PP-1400	----
Mint	#PP-1100	#PP-1500	----
Cherry	#PP-1200	#PP-1600	----
Bubble Gum	----	----	#PP-1800

\*Assorted contains: cherry, mint, grape and raspberry  
(assortment could vary)

## Defend Plus Prophy Angles



- Contra design reduces hand fatigue, micro traumas and allows for greater access
- Disposable • Latex Free
- Smooth running, vibration free
- Flared cup design reduces splatter

100 per box

#PA-6000 - firm (Purple)

#PA-7000 - soft (Green)



Lowest Priced  
Premium Contra Angle  
on the Market!



## Pit & Fissure Sealant - Light Curing

Covers and seals pits and fissures in occlusal surfaces of posterior teeth to aid resistance to caries on those surfaces. Its very low viscosity resin system has a white opaque color to help identify retention at recall appointments. The direct delivery syringe system provides a quick and easy application method. The resin cures quickly to form a hard and retentive sealant.

Contains: 4 - 1GM Syringes LC Opaque Sealant, 1 - 2ml Syringe  
Etchant Gel, 6 - Black/Blue Metal Pre-Bent Syringe Tips (12 total)

#PF-1000



## Etch Gel - 38%

A 38% phosphoric acid gel for etching tooth enamel. Blue color is easily distinguishable from the tooth structure. Pre-bent tips allow for precise placement.

Contains: 2 - 2ml Syringes Etchant Gel,  
20 - Blue Metal Pre-Bent Syringe Tips

#PF-2000



## Defend Toothbrushes

- Same quality of major brands at a fraction of the cost
- Individually wrapped
- Assorted colors



#TB-2500

#TB-2000

#TB-1500

#TB-1000

72 per box/6 boxes per case


#TB-1000 Premium Adult Toothbrush

#TB-1500 Premium Adult Toothbrush with rubber grip

#TB-2000 Youth Toothbrush

#TB-2500 Children's Toothbrush

# Disposables




## Disposable Aseptic Bib Holders



Disposable (Paper) Bib Holders. No Sharp Edges or Cold Chain. Stretch to fit, Self-Adhesive ends.

250 per box  
#PB-8500



## Disposable 4-Ply Patient Towels



DEFEND's new 4-ply patient towels are impervious to moisture and saliva. Constructed with a new 3-ply embossed tissue and 1-ply poly, these towels provide maximum absorption and patient protection.

13"x18" - 500 per case  
#PB-9001 Blue



## Disposable Patient Towels



The unique technique used to bond the layers of tissue and poly securely together eliminates separation of the layers, and reinforced edges provide added strength and durability. 2-ply tissue and 1-ply poly for absorption and patient protection

13"x18" - 500 per case  
#PB-8000 White  
#PB-8001 Blue  
#PB-8002 Mauve  
#PB-8003 Green  
#PB-8004 Lavender  
#PB-8005 Yellow



## Disposable Paper Tray Covers



Made of durable construction, these paper tray covers provide low absorbency and fit trays just right. They feature rounded corners and come in the most popular "B" size - 8½" x 12¼" (21.6 cm x 31.1 cm).

1000 per case  
#TC-9000 White  
#TC-9001 Blue  
#TC-9002 Mauve  
#TC-9003 Green  
#TC-9004 Lavender  
#TC-9005 Yellow



## 6" Adult - Disposable Tongue Depressors

These 100% wood tongue depressors are smooth and splinter-free. Non-sterile.

500 per box  
#TD-9000



## Green Seal® Certified Paper Products



- 100% recycled paper products
- Chlorine-free, with no added pigments, dyes or fragrances
- 2-ply facial and toilet tissue are soft & strong
- C-fold towels are soft and ultra absorbent

**C-Fold Towel**  
#CF-1000 - 10 1/8" x 13"  
200 sheets/pkg - 12 per case

**Facial Tissue**  
#FT-1000 - 8" x 8" - 100 count  
30 per case

**Toilet Tissue**  
#BT-1000 - 500 sheets per roll  
48 per case




## Disposable Drinking Cups



These disposable drinking cups have a special embossed grip that provides greater strength and makes dispensing easier. They have a lustrous high-gloss finish and are available in five designer colors. Recyclable.

5oz/148 ml - 1000 per case  
#DC-7000 White  
#DC-7001 Blue  
#DC-7002 Mauve  
#DC-7003 Green  
#DC-7004 Lavender  
#DC-7007 Yellow



## Cotton Tipped Applicators

Non-sterile, made from all wood dowels with very tight 100% cotton bud to minimize linting and prevent shredding. Comes sealed in autoclavable bag with indicator.

1000 per box  
#CS-3000 3" dowel  
#CS-6000 6" dowel



## Cotton Rolls



High quality, soft, made from 100% virgin cotton, are non-linting, non-sterile, pliable and 20% more absorbent than the leading brand. They hold together better, maintain their shape when wet and will not adhere to mucous membranes.


#2, medium 1.5" x .375" (3/8")  
2000 per box  
#CS-0200



## Headrest Covers

Made with a soft, comfortable textured tissue outer covering and a poly backing to help prevent moisture soak-through. Fluid-resistant and disposable.

White - 500 per case  
#HC-2001 10" x 10"  
#HC-1001 10" x 13"



## 2" x 2" Cotton Filled Gauze Sponges

These cotton filled gauze sponges contain 20% more cotton and are more absorbent. These pure white, 8-ply sponges are shred-resistant and made of 100% cotton filled USP gauze.

5000 per case  
#CS-0100 Non-Sterile  
#CS-0101 Sterile



## Non-Woven Sponges

A great alternative to cotton filled sponges, these are made from a 4-ply, 32 gram rayon/poly blend. Ideal for wound cleansing, and as secondary dressings for added absorbency. Non-linting/non-sterile.

#NW-0200 2" x 2" - 5000/cs  
#NW-0300 3" x 3" - 4000/cs  
#NW-0400 4" x 4" - 2000/cs





## Multi-Purpose Disposable Dappen Dishes



Designed to prevent cross contamination for a variety of procedures. A quick and convenient vessel for prophylactic paste, composite materials, acrylics, cements and bonding agents. Each dish is made of high quality plastic that is impervious to most dental materials. Assorted colors.

1000 per box  
#DD-5000



## DEFEND+PLUS Pre-Bent Dispensing Tips

Pre-bent to the proper angle for placing etchants, sealants, cements and other flowable materials. German made & final step washed and blown dry to avoid sediment within the metal tubes. Color coded represent gauges.

#PN-1227 Pink [Liners & Bases (18 ga.)]  
#PN-1228 Yellow [Sealant & Composites - 20ga.]  
#PN-1229 Black [Composites & Cements - 19ga.]  
#PN-1230 Grey [Sealant & Cement - 22ga.]  
#PN-1231 Blue [Etchants - 25ga.]

## Defend 20% Benzocaine Anesthetic Gels



- Fast acting with no systemic absorption
- Provides temporary relief of pain during dental procedures
- Comforts patients with minor mouth irritations
- No bitter after taste

1oz. (30g)

#TA-5000 - Strawberry #TA-5005 - Grape  
#TA-5001 - Cherry #TA-5006 - Mango  
#TA-5002 - Pina Colada #TA-5007 - Raspberry  
#TA-5003 - Mint #TA-5008 - Banana  
#TA-5004 - Bubble Gum

## Plastic Hub/Self-Threading Dental Needles



- Needle made of stainless steel; polypropylene hub with aluminum insert
- Disposable, sterile needles for use on standard 1.8mL dental syringes
- Triple-beveled needle features our unique arrow indicator to help orient the bevel position during injections
- Siliconized to help insertion of needle
- Pre-threaded plastic hub for simple, straight attachment to the syringe
- Color coding for easy needle diameter identification


Yellow #DN-2000 27L #DN-2100 27S  
Blue #DN-3000 30S #DN-3100 30XS



## Bendable - Tapered Applicator Brushes

Makes it easier to apply whitening agents, etchants, bonding agents, sealants, desensitizers, adhesives and resins. The brush's chemical resistant fibers hold up to the harshest dental materials and they allow for greater carrying capacity and better smoothing and spreading.

100 per tube  
#BB-1440 Blue  
#BB-1450 Pink



## Disposable - 7 Hole Air/Water 3-Way Syringe Tips

Disposable tip reduces the potential for cross-contamination via air/water syringe. Delivers precise separation of air and water for consistent dry air, and is ideal for sealant applications and restorative procedures.

- CLEAR disposable tip • Consistent dry air
- Smooth exterior and polished ends
- Rigid – provides secure retraction • Latex free

250 per bag  
#AW-6000

Sani-Tip® is a registered trademark of DENTSPLY International Inc.



## Disposable Metal Core Air/Water 3-Way Syringe Tips

- Replaces existing metal (autoclavable) tips
- Inner core water chamber
- Secure o-ring groove eliminates leakage
- Strong enough to allow for retraction
- Unique quad core design ensures dry air
- Prevents cross contamination
- Delivers superior air/water spray

150 per bag  
#AW-4000

## Plastic Air/Water 3-Way Syringe Tips



Excellent for the replacement of existing metal tips and has the identical angle. The plastic material designed tip won't damage the O-ring like metal tips can. An O-ring groove secures tip for patient safety and eliminates leakage. Unique quad-core design ensures instant dry air and prevents cross contamination.

250 per bag 1500 per bag  
#AW-1000 - Clear #AW-3500 - Blue  
#AW-2000 - White  
#AW-3000 - Blue

## Disposable • 12cc Curved Utility Syringe



- Multipurpose • Non-Sterile
- Use for Impression Materials and Irrigation

50 per box  
#IS-1200



# Barrier Products



## Barrier Film

Universal infection control barrier with a low tack adhesive backing with "non-stick" edge. It is perfect for use in areas that are difficult to clean, disinfect or sterilize. Perforated 4" x 6" sheets leave behind no residue.

1200 sheets per roll

In Dispenser Box

#BF-1500 Clear


#BF-2500 Blue

Rolls Only

#BF-1600 Clear

#BF-2600 Blue

## Barrier Film Dispenser



(barrier film not included)

This reusable dispenser reduces waste of dispenser boxes while minimizing office clutter. Use with item's #BF-1500/1600 and #BF-2500/2600.

1 Unit

#BF-1800



## iPad/Tablet Sleeves

iPads and tablets have become more prevalent for staff use and patient communication. These invaluable tools are also the most subject to cross-contamination. These barriers are specifically designed to protect touch screen products while not interfering with their performance. Protect yourself and these valuable technology items from splash, splatter, over sprays and other risk.



NEW!

100 per box

#BF-7800



## Syringe Sleeves

Fits most 3-way air/water syringes and has a pre-punched opening making application easy. They can also be used to cover handpieces, saliva ejectors, ultrasonic scalers and more.



500 per box

#BF-3000 - 2½" x 10"



## Handpiece Sleeves

Made to fit over most handpiece instruments. Made of high quality clear plastic.



500 per box

#IN-7006 - 1½" x 8"



## Fitted Digital Sensor Sleeves

- Precise, snug fit to sensor size resulting in no additional plastic in patient's mouth
- Twice the thickness of conventional sleeves
- Reduced tearing • Soft and rounded edges
- Latex-free

500 per box

Size 1: 1.125"(W) x 9.25"(L)

x 1.25" Head Size)

Fits: Kodak, Schick, Gendex, EVA & Planmeca

#BF-8200

Size 2: Fits: Kodak 6100, Schick, Gendex GX-S, EVA and Planmeca

#BF-8250

Snug, tight fit to sensor.



## Chair Sleeves

Our Chair Sleeves are made of a clear, highly durable plastic and are designed to easily fit over most chairs.

225 per box

1/2 Chair - #BF-4000 - 27½" x 24"

125 per box

Full Chair - #BF-4500 - 29" x 79"



## Headrest Covers

Designed to fit over the top of headrests. Made of clear plastic and are easy to apply.

250 per box

#BF-9000 - 11" x 9½" x 2"

#BF-9500 - 14" x 9½" x 2"



## X-Ray Head Sleeves

Our X-Ray Head Sleeves are designed to fit over most X-Ray heads. The barrier is made of clear plastic.



250 per box

#BF-6000 - 15" x 26"



## Tray Sleeves

Tray Sleeves cover the entire instrument tray and are made of clear plastic.



500 per box

#BF-7000 - 10½" x 14"

#BF-7500 - 11½" x 16"



## Digital X-Ray Sensor Sleeves

Protect your investment in digital imaging with extra soft, tear-resistant, leak-proof construction.

500 per box

#BF-8000

1½" x 8"

Trophy, Gendex & Diagnodent

#BF-8100

2½" x 10"

Various systems

#BF-8500

1½" x 8"

Trophy & Gendex



NEW!



## Defend Plus Glider Barrier Envelopes

The perfect barrier envelope to protect your phosphor plates from cross contamination. Features a unique "GLIDER" loader/sealer allowing you to quickly and easily "glide" your PSP's into the envelope. Upon peeling off the loader's tab, the envelope practically seals itself. When the x-rays are completed, the center seam allows you to tear open the envelope and start the developing process without ever touching the plate. Can be used with virtually any system on the market.

100 per box

#BF-8750 - #0

100 per box

#BF-8850 - #1

300 per box

#BF-8950 - #2



## Barrier Envelopes for phosphor storage plate



Protect your investment in digital imaging. Features a tear-away thumb-notch opening that makes it easy to remove plates for processing. Made of soft, supple vinyl with rounded corners to maximize patient comfort.

#BF-8700 #0 100/pkg

#BF-8800 #1 100/pkg

#BF-8900 #2 300/pkg



## Light Handle Sleeves

Light Handle Sleeves are made of clear plastic and are made to fit most light handles.

500 per box

#BF-5000 - 4" x 5¾"





**New** sure grip is made of autoclavable silicone.



## MirrorLite™ Illuminated Mouth Mirror

The portable MirrorLite™ offers superb, shadow and glare-free illumination of the oral cavity and tooth structures, without cables or the need to readjust overhead lights! Perfect for a general exam as well as the most detailed hygiene procedure.

- Uses pure white LED light
- Titanium coating resists corrosion and tarnishing
- Knurled finish provides a better grip
- Ergonomic, compact design reduces eye and back strain

### #IN-7003 - Introductory Kit

(contains: 1 MirrorLite, 50 Disposable Mirrors, 6 batteries and 50 Barrier Sleeves)

### #IN-7004 - #5 Refill Mirrors

(contains: 50 Mirrors, 2 batteries and 50 Barrier Sleeves)

### #IN-7005 - MirrorLite only

### #IN-7006 - Barrier Sleeves - 500 per box

### #IN-7007 - Batteries - 3/pk

Introductory Kit - #IN-7003



## Micro Torch & Butane Refill

- Lightweight & compact
- Instant electronic ignition system
- Refillable butane gas chamber
- Includes a removable base for free-standing or handheld use

#MT-9000 Torch

#MT-9001-Butane Refill

5.8 oz. (16.5g)



## Handpiece Cleaner & Lubricant

- Synthetic hydrocarbons
- Increases range of high 485° F to -85° F
- Metered, non-aerosol spray
- Prevents O-rings from drying out
- Low evaporation
- Prevents oxidation & corrosion
- Non-toxic and non-flammable
- Extends the life of any handpiece

2 x (2 oz.) bottles per kit

#SO-9700



## Restorative



## Temporary Crown & Bridge Material

DEFEND Temporary Crown & Bridge Material is a syringeable bis-acryl composite for chairside provisional restorations. DEFEND Temporary Crown & Bridge Material is fast and accurate, especially when used in conjunction with a quality vinyl polysiloxane impression material such as DEFEND Super Hydrophilic VPS (Monophase). DEFEND Temporary Crown & Bridge Material is available in 6 shades: A1, A2, A3, A3.5, B1, and Bleach BL. All are delivered in 76g (50ml) cartridges designed to fit on a 10:1 style automix gun.

Contains: 76gm 10:1 cartridge & 10 T-Mixer tips (material saving tip)

#CB-9000 - A1

#CB-9001 - A2

#CB-9005 - A3

#CB-9002 - A3.5

#CB-9003 - B1

#CB-9004 - Bleach



## Disposable Jackets & Lab Coats



**The Latest Technology in Protection & Comfort!**

*Provides 28% better filtration efficiency compared to leading competitor.\**

- Waffle Patterned Fabric
- Soft 3-Layer "SMS" Fabrics
- Fluid Resistant, Breathable
- Snap Front, Knit Cuffs & Collar
- 1 Breast Pocket & 2 Hip Pockets
- Autoclavable\*\* & Washable
- Lightweight, Latex Free
- Static Free - No Clinging



Ceil Blue



Light Pink



Full Length Lab Coat

### Hip Length Jacket - 10/Bag

Size	Ceil Blue	Sky Blue	Purple	White
Small	#SJ-3630CBS	#SJ-3630SBS	#SJ-3630PPS	#SJ-3630WS
Medium	#SJ-3630CBM	#SJ-3630SBM	#SJ-3630PPM	#SJ-3630WM
Large	#SJ-3630CBL	#SJ-3630SBL	#SJ-3630PPL	#SJ-3630WL
XLarge	#SJ-3630CBXL	#SJ-3630SBXL	#SJ-3630PPXL	#SJ-3630WXL
XXLarge	#SJ-3630CBXXL	-	-	-

### Full Length Lab Coat - 10/Bag

Size	Ceil Blue	Sky Blue	Purple	White
Small	#SG-3660CBS	#SG-3660SBS	#SG-3660PPS	#SG-3660WS
Medium	#SG-3660CBM	#SG-3660SBM	#SG-3660PPM	#SG-3660WM
Large	#SG-3660CBL	#SG-3660SBL	#SG-3660PPL	#SG-3660WL
XLarge	#SG-3660CBXL	#SG-3660SBXL	#SG-3660PPXL	#SG-3660WXL
XXLarge	#SG-3660CBXXL	-	-	-

**NEW!**

### FITTED - 10/Bag

#### Hip Length Jacket

#### Full Length Lab Coat

Size	Light Pink	Light Pink
Small	#SG-3630LPS	#SG-3660LPS
Medium	#SG-3630LPM	#SG-3660LPM
Large	#SG-3630LPL	#SG-3660LPL
XLarge	#SG-3630LPXL	#SG-3660LPXL

Close up of SMS waffle pattern fabric.

**NOW AVAILABLE IN WHITE!**

THE DENTAL ADVISOR

++++ 1/2

#### Available Colors:

Ceil Blue	Sky Blue	Purple	White	Light Pink

\*According to Nelson Laboratories, a leading USA medical testing facility 12/14/10  
 \*\*Note: close snaps when autoclaved

## Disposable Lab Coat



DEFEND has improved the quality of its lab coats to make them more comfortable for healthcare professionals. The new coats are more breathable, while retaining the same fluid-resistance they've always had. These disposable/reusable lab coats feature a front pocket, easy-to-use snap closures, knit collar and cuffs. The unisex, knee-length blue coats are available in 6 sizes.

*10 per bag*

*#SG-9003 Small*

*#SG-9004 Medium*

*#SG-9005 Large*

*#SG-9006 X-Large*

*#SG-9007 XX-Large*

*#SG-9008 XXX-Large*



## Disposable Tie-Back Protective Gown



These full-length gowns are made from a fluid resistant, cool and comfortable material to provide both protection and comfort. They are an inexpensive alternative to lab coats, yet have key features such as waist ties and elastic cuffs. The unisex, knee-length blue coats are available in 2 sizes.

*10 per bag*

*#SG-1000 Medium (40" length)*

*#SG-1001 Large (48" length)*



## Clear Protective Eyewear



This economical eyewear offers excellent protection in the operator or lab. The one piece lens is made of a high impact, scratch resistant clear polycarbonate material. They provide non-glare and optically clear protection.

*1 per box*

*#GL-2021*



## U.V. Protective Eyewear



This eyewear filters out 98% of dangerous blue light frequencies and has total protection against UV light. The one piece lens is made of a high impact, scratch resistant polycarbonate material.

*1 per box*

*#GL-2020*



## ORTHO DENTAL BOXES



### Retainer Boxes

Crush proof retainer boxes have a solid locking mechanism and easily cleans with soap and water.

*3" x 2½" x 1" deep*

*12/box in 4 assorted colors (blue, white, red and Yellow)*

*#OB-2000*



### Denture Boxes

High-impact plastic, 1 piece construction and contoured corners for ease of cleaning. Tight fitting lid holds 1 or 2 dentures.

*3" x 2½" x 2" deep*

*12/box in 4 assorted colors (blue, white, mauve and beige)*

*#OB-1000*

## Mydent Affiliations:



Mydent International is a proud supporter of OSAP, the Organization for Safety, Asepsis and Prevention. OSAP advocates for the safe and infection-free delivery of oral care. OSAP targets its message to dental and other healthcare professionals, academia, consultants, researchers and organizations that include non-governmental organizations (NGOs), government agencies, and even manufacturing and distribution companies -- it's a trusted resource for everyone who is responsible for safe dental care...because safety matters! Learn more at [www.osap.org](http://www.osap.org).



Mydent International is an accredited business with the Better Business Bureau (BBB) and has a prestigious A+ rating. Since 1985, we have committed ourselves to providing ethical business practices, honest advertising and safeguarding customer privacy.



Mydent International is proud to be a longtime member of the Dental Trade Alliance. The DTA is an association of companies that provide dental equipment, supplies, materials and services to dentists and other oral care professionals. Our member companies are distributors, dental laboratories & manufacturers located in the United States, Canada and Mexico. By providing the best equipment, materials and services to dentists and oral care professionals, we are partners in improving the oral health of everyone.



With Autism rates estimated as high as 1 in 91 children, Mydent feels compelled to take action. That is why a portion of the profits from every DEFEND product sold will go to Autism Speaks. Founded in 2005, Autism Speaks has quickly grown to become the nation's largest and most effective autism science and advocacy organization. They are dedicated to funding global biomedical research into the causes, prevention, treatments, and cure for autism. Mydent is proud to play a small part in helping them to change the future for all who struggle with autism spectrum disorders. We urge you to learn more at [www.autismspeaks.org](http://www.autismspeaks.org).



Mydent is also becoming an industry leader in environmentally responsible manufacturing and distribution practices. We are committed to a reduction of our carbon footprint, using earth-friendly materials and processes wherever possible. From recyclable and reduced packaging, high order practices, reducing emissions as well as "green" practices at our own headquarters, we are committed to a philosophy of being environmentally responsible...Today and Tomorrow.

### Our Brand Promise

**"To provide the healthcare professional with the highest quality infection control products, disposables, preventatives and impression material systems at affordable prices, supported by Superior Service and 100% Customer Satisfaction."**

# Mydent International

80 Suffolk Court, Hauppauge, NY 11788  
800.275.0020 • FAX: 631.434.7750  
Ph: 631.434.3190 • EMAIL: [sales@defend.com](mailto:sales@defend.com)



[www.defend.com](http://www.defend.com)

## JOIN OSAP.ASAP

The Organization for Safety, Asepsis and Prevention (OSAP) is the world's leading advocate and trusted advisor for the safe and infection-free delivery of oral health care.

The Safest Dental Visit™ is OSAP's latest infection control educational program. As an authority on infection control, OSAP has developed this program to help dental practices:

- Offer patients The Safest Dental Visit™ - elevating their confidence in the practice
- Engage staff in a "culture of safety"
- Enhance the overall image of the practice

Mydent  
International

IS A PROUD MEMBER OF



It's Like Having an Infection Control Partner! Visit: [www.osap.org](http://www.osap.org)